

Overview wet cart

Restart button



FACS Flow

- FACSFlow is supplied in 20 l cartons.
- New cartons can be found in WAD-K-126.

Waste

- Waste is collected in 10 l plastic tanks.
- An empty waste tank should be next to the wet cart. If not check the sink and empty the tank with decontaminated waste into the sink.

Exchange FACSFlow or Waste

FACSFlow

- The “1.5 liters” lamp indicates low sheath fluid levels; no user intervention is required.
- When the alarm sounds and the “Empty” lamp is on, following user action is required.
 - Mute the alarm sound by pressing “alarm” on the wet cart.
 - Finish measuring your current sample tube and remove it from the SIT arm.
 - Switch the instrument to “Standby” leaving the sample arm open.
 - Replace FACSFlow carton with full one.
 - Place the empty carton in the hallway next to class door at the elevator.
 - Press the “**Restart button**” on the fluidics cart for >1 sec and continue working (all red lights should turn off).

Waste Full

- When the alarm sounds and the “Waste Full” lamp is on, following user action required.
 - Mute the alarm sound by pressing “alarm” on the wet cart.
 - Finish measuring your current sample tube and remove it from the SIT arm.
 - Switch the instrument to “Standby” leaving the sample arm open.
 - Replace the full waste tank with the empty one.
 - Bring the full waste tank to sink and decontaminate it by adding 100 ml Javel (blue canister in the cupboard below the sink).
 - Press the “**Restart button**” on the fluidics cart and continue working.